

The Implant Essential



Clinically Proven
Up To **145%**
More Effective for
Implant Patients
vs. String Floss.¹



Clinically Proven Defense Against Peri-Implantitis

Implants are today's treatment of choice for missing dentition, and a hygiene challenge. Up to half of implant patients now have peri-implantitis, which may result in implant failure.²

Given the depth of most implant pockets, traditional brushing and flossing is insufficient. That's why we developed the patented Plaque Seeker® Tip. It allows 360-degree cleaning from the gingival margin to the base of the pocket. And new clinical research from Tufts University shows the Waterpik® Water Flosser with the Plaque Seeker® Tip is up to 145% more effective than traditional string floss for improving gingival health around implants.¹

So help your patients protect their investment with the Waterpik® Water Flosser - the essential tool that no implant patient should be without.

"For all implant restorations, I recommend the Waterpik® Water Flosser. Its uniquely designed tips are perfect for delivering irrigation around the challenging contours of an implant restoration."

Chris Salierno, DDS

To see Clinical Research Abstracts and
Published Studies go to: Professional.Waterpik.com

waterpik®

WATER FLOSSER

The Easy and More
Effective Way to Floss™

1 - Clinical Study, Tufts University, 2013. See waterpik.com for details.
2 - J. Clin Periodontol 2008; 35 (Suppl. 8): 282-285. ©2013 Water Pik, Inc.